

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-African clawed frog Nup54 NANOBODY [®] reg, protein (SAA0873)
PRODUCT REFERENCE	ARO-A16282
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	Nup54

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Xenopus laevis</i> (African clawed frog)
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Nup54, anti-nucleoporin 54kDa S homeolog isoform X2, anti-nup54.S, anti-nup54