

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human GPHN/Gephyrin NANOBODY®, protein (SAA1333)
PRODUCT REFERENCE	ARO-A16284
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	GPHN/Gephyrin

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-MPT Mo-transferase, anti-GPHN, anti-GPH, anti-Gephyrin, anti-KIAA1385, anti-Domain G, anti-Domain E, anti-MPT adenylyltransferase