



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-SARS-CoV S/Spine glycoprotein NANOBODY®, protein (SAA1322)
PRODUCT REFERENCE	ARO-A16289
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	S/Spine glycoprotein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Severe acute respiratory syndrome coronavirus (SARS-CoV)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Spine glycoprotein, anti-S glycoprotein, anti-E2, anti-Peplomer protein, anti-Spine protein S1, anti-Spine protein S2, anti-Spine protein S2', anti-S, anti-2, anti-Severe acute respiratory syndrome coronavirus (SARS-CoV)

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY