

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mouse H2-Ab1/MHCII NANOBODY [®] protein (SAA1323)
PRODUCT REFERENCE	ARO-A16291
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	H2-Ab1/MHCII

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-H-2 class II histocompatibility antigen, anti-A beta chain, anti-H2-Ab1, anti-H2-iabeta