

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-SARS-CoV-2 NP/Nucleoprotein NANOBODY [®] , protein (SAA0960)
PRODUCT REFERENCE	ARO-A16294
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	NP/Nucleoprotein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Nucleoprotein, anti-N, anti-Nucleocapsid protein, anti-NC, anti-Protein N, anti-N, anti-Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)