

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Bacillus anthracis PA/pagA/PA-83 NANOBODY <sup>®</sup> , protein (SAA1225)
<b>PRODUCT REFERENCE</b>	ARO-A16302
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VHH-8His-Cys-tag
<b>TARGET</b>	PA/pagA/PA-83

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Bacillus anthracis</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	anti-Protective antigen, anti-PA, anti-Anthrax toxins translocating protein, anti-PA-83, anti-PA83, anti-Protective antigen PA-20, anti-PA-20, anti-PA20, anti-Protective antigen PA-63, anti-PA-63, anti-PA63, anti-pagA, anti-pag, anti-pXO1-110, anti-BXA016