

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HCMV/HHV5 HHRF3/US28 NANOBODY [®] , protein (SAA0899)
PRODUCT REFERENCE	ARO-A16307
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	HHRF3/US28

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human cytomegalovirus (strain Towne) (HHV-5) (Human herpesvirus 5)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-G-protein coupled receptor homolog US28, anti-HHRF3, anti-US28, anti-Human cytomegalovirus (strain Towne) (HHV-5) (Human herpesvirus 5)