

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Clostridium botulinum botA/BOTOX NANOBODY <sup>®</sup> , protein (SAA0923)
<b>PRODUCT REFERENCE</b>	ARO-A16308
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VHH-8His-Cys-tag
<b>TARGET</b>	botA/BOTOX

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Clostridium botulinum</i> (strain Hall / ATCC 3502 / NCTC 13319 / Type A)
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	anti-Botulinum neurotoxin type A, anti-BoNT/A, anti-Bontoxilysin-A, anti-BOTOX, anti-Botulinum neurotoxin A light chain, anti-LC, anti-3.4.24.69, anti-Botulinum neurotoxin A heavy chain, anti-HC, anti-botA, anti-bna, anti-CBO0806, anti-CLC_0862