

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HRSV F/Fusion glycoprotein F0 NANOBODY [®] , protein (SAA1037)
PRODUCT REFERENCE	ARO-A16309
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	F/Fusion glycoprotein F0

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human respiratory syncytial virus A (strain A2)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Fusion glycoprotein F0, anti-Fusion glycoprotein F2, anti-F2, anti-p27, anti-Intervening segment, anti-Pep27, anti-Peptide 27, anti-Fusion glycoprotein F1, anti-F1, anti-F, anti-Human respiratory syncytial virus A (strain A2)