

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human BACE2 NANOBODY [®] , protein (SAA1031)
PRODUCT REFERENCE	ARO-A16316
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	BACE2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Aspartyl protease 1, anti-Membrane-associated aspartic protease 1, anti-Beta-secretase 2, anti-ASP21, anti-BACE2, anti-Asp 1, anti-ALP56, anti-Beta-site APP cleaving enzyme 2, anti-Aspartic-like protease 56 kDa, anti-Beta-site amyloid precursor prote