

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human IL17RB NANOBODY [®] , protein (SAA1267)
PRODUCT REFERENCE	ARO-A16322
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	IL17RB

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-IL-17B receptor, anti-CRL4, anti-IL17RB, anti-Cytokine receptor-like 4, anti-IL-17RB, anti-IL17BR, anti-IL-17 receptor B, anti-Interleukin-17B receptor, anti-EVI27, anti- Interleukin-17 receptor B, anti-IL17Rh1, anti-IL-17Rh1, anti-IL-17 receptor homo