

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human FOLH1/PSMA NANOBODY [®] , protein (SAA1352)
PRODUCT REFERENCE	ARO-A16323
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	FOLH1/PSMA

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-NAALAD1, anti-Prostate-specific membrane antigen, anti-NAALADase I, anti-PSMA, anti-PSM, anti-Pteroylpoly-gamma-glutamate carboxypeptidase, anti-Glutamate carboxypeptidase 2, anti-N-acetylated-alpha-linked acidic dipeptidase I, anti-Glutamate carboxy