

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD47/MER6 NANOBODY [®] , protein (SAA1315)
PRODUCT REFERENCE	ARO-A16329
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD47/MER6

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Leukocyte surface antigen CD47, anti-CD47, anti-Antigenic surface determinant protein OA3, anti-IAP, anti-Integrin-associated protein, anti-Protein MER6, anti-MER6