

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT NAME | Anti-Human Adiponectin/ACRP30/ADIPOQ NANOBODY [®] , protein (SAA1153) |
| PRODUCT REFERENCE | ARO-A16330 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | VHH-8His-Cys-tag |
| TARGET | Adiponectin/ACRP30/ADIPOQ |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Homo sapiens (Human)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4) |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | anti-Adiponectin, anti-ACRP30, anti-GBP28, anti-Adipose most abundant gene transcript 1 protein, anti-ACDC, anti-APM1, anti-ADIPOQ, anti-apM-1, anti-30 kDa adipocyte complement-related protein, anti-Adipocyte complement-related 30 kDa protein, anti-Adipoc |