

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human ASC/TMS1/PYCARD NANOBODY&reg, protein (SAA1187)
<b>PRODUCT REFERENCE</b>	ARO-A16331
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VHH-8His-Cys-tag
<b>TARGET</b>	ASC/TMS1/PYCARD

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	anti-TMS1, anti-CARD5, anti-Caspase recruitment domain-containing protein 5, anti-Apoptosis-associated speck-like protein containing a CARD, anti-PYD and CARD domain-containing protein, anti-ASC, anti-PYCARD, anti-hASC, anti-Target of methylation-induced