

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ASC/TMS1/PYCARD NANOBODY®, protein (SAA1025)
PRODUCT REFERENCE	ARO-A16334
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	ASC/TMS1/PYCARD

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-TMS1, anti-CARD5, anti-Caspase recruitment domain-containing protein 5, anti-Apoptosis-associated speck-like protein containing a CARD, anti-PYD and CARD domain-containing protein, anti-ASC, anti-PYCARD, anti-hASC, anti-Target of methylation-induced