

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human TRIM28 NANOBODY [®] reg, protein (SAA1235)
PRODUCT REFERENCE	ARO-A16335
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	TRIM28

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-TIF1-beta, anti-Tripartite motif-containing protein 28, anti-KRAB-associated protein 1, anti-TRIM28, anti-RING-type E3 ubiquitin transferase TIF1-beta, anti-KRAB-interacting protein 1, anti-KRIP-1, anti-E3 SUMO-protein ligase TRIM28, anti-RNF96, anti