

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ADAMTS5 NANOBODY [®] reg, protein (SAA1219)
PRODUCT REFERENCE	ARO-A16340
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	ADAMTS5

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-ADAMTS11, anti-ADAM-TS 11, anti-A disintegrin and metalloproteinase with thrombospondin motifs 11, anti-ADAMTS-11, anti-ADMP2, anti-ADAM-TS 5, anti-ADAMTS-5, anti-ADMP-2, anti-ADAMTS5, anti-Aggrecanase-2, anti-ADAM-TS5, anti-A disintegrin and metallo