

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ADAM9 NANOBODY [®] , protein (SAA1216)
PRODUCT REFERENCE	ARO-A16341
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	ADAM9

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Disintegrin and metalloproteinase domain-containing protein 9, anti-Meltrin-gamma, anti-ADAM9, anti-MDC9, anti-MLTNG, anti-ADAM 9, anti-Metalloprotease/disintegrin/cysteine-rich protein 9, anti-Cellular disintegrin-related protein, anti-KIAA0021, ant