

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human BAIAP2 NANOBODY [®] reg, protein (SAA1334)
PRODUCT REFERENCE	ARO-A16346
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	BAIAP2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-BAI-associated protein 2, anti-Insulin receptor substrate p53, anti-IRS-58, anti-BAI1-associated protein 2, anti-IRSp53, anti-Insulin receptor substrate p53/p58, anti-Fas ligand-associated factor 3, anti-BAIAP2, anti-FLAF3, anti-Insulin receptor subs