

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human WFDC2/WAP5 NANOBODY®, protein (SAA1359)
PRODUCT REFERENCE	ARO-A16364
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	WFDC2/WAP5

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Epididymal secretory protein E4, anti-Major epididymis-specific protein E4, anti-HE4, anti-WFDC2, anti-WAP5, anti-WAP four-disulfide core domain protein 2, anti-Putative protease inhibitor WAP5