

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |  |
|--------------------------|--|
| <b>PRODUCT NAME</b>      | Anti-Human SLC10A1 NANOBODY <sup>®</sup> protein (SAA1145) |
| <b>PRODUCT REFERENCE</b> | ARO-A16368   |
| <b>CLONALITY</b>         | Monoclonal Antibody  |
| <b>ISOTYPE</b>           | VHH-8His-Cys-tag   |
| <b>TARGET</b>            | SLC10A1  |

### PRODUCT INFORMATION

---

|                                |   |
|--------------------------------|---|
| <b>ORIGIN SPECIES</b>          | <i>Homo sapiens (Human)</i>   |
| <b>FORM</b>                    | Liquid  |
| <b>BUFFER</b>                  | 0.01M PBS (pH 7.4)  |
| <b>STABILITY &amp; STORAGE</b> | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| <b>ALIAS</b>                   | anti-Cell growth-inhibiting gene 29 protein, anti-Sodium/taurocholate cotransporting polypeptide, anti-NTCP, anti-SLC10A1, anti-Sodium/bile acid cotransporter, anti-Na(+)/taurocholate transport protein, anti-Solute carrier family 10 member 1, anti-Na(+)/b   |