

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD202b/TEK/TIE-2 NANOBODY [®] protein (SAA1254)
PRODUCT REFERENCE	ARO-A16370
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD202b/TEK/TIE-2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-VMCM, anti-TIE2, anti-Tunica interna endothelial cell kinase, anti-CD202b, anti-Tyrosine-protein kinase receptor TIE-2, anti-TEK, anti-Angiopoietin-1 receptor, anti-VMCM1, anti-Endothelial tyrosine kinase, anti-Tyrosine-protein kinase receptor TEK, a