

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human SLC10A1 NANOBODY [®] protein (SAA1065)
PRODUCT REFERENCE	ARO-A16371
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	SLC10A1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Cell growth-inhibiting gene 29 protein, anti-Sodium/taurocholate cotransporting polypeptide, anti-NTCP, anti-SLC10A1, anti-Sodium/bile acid cotransporter, anti-Na(+)/taurocholate transport protein, anti-Solute carrier family 10 member 1, anti-Na(+)/b