

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human AKT1 NANOBODY [®] reg, protein (SAA1231)
PRODUCT REFERENCE	ARO-A16376
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	AKT1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-RAC-alpha serine/threonine-protein kinase, anti-RAC, anti-RAC-PK-alpha, anti-PKB alpha, anti-Proto-oncogene c-Akt, anti-Protein kinase B alpha, anti-Protein kinase B, anti-AKT1, anti-PKB