

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human BACE1 NANOBODY [®] , protein (SAA1336)
PRODUCT REFERENCE	ARO-A16381
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	BACE1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-BACE1, anti-Memapsin-2, anti-KIAA1149, anti-Asp 2, anti-Beta-site amyloid precursor protein cleaving enzyme 1, anti-BACE, anti-Beta-site APP cleaving enzyme 1, anti-Beta-secretase 1, anti-Aspartyl protease 2, anti-ASP2, anti-Membrane-associated aspar