

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD197/CCR7 NANOBODY®, protein (SAA1349)
PRODUCT REFERENCE	ARO-A16382
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD197/CCR7

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-CDw197, anti-EVI1, anti-CMKBR7, anti-CCR7, anti-C-C CKR-7, anti-CD197, anti-CCR-7, anti-Epstein-Barr virus-induced G-protein coupled receptor 1, anti-MIP-3 beta receptor, anti-EBI1, anti-CC-CKR-7, anti-BLR2, anti-EBV-induced G-protein coupled recepto