

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CASP3/Caspase-3 NANOBODY [®] reg, protein (SAA1223)
PRODUCT REFERENCE	ARO-A16383
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CASP3/Caspase-3

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Caspase-3, anti-CPP-32, anti-CPP32, anti-CASP-3, anti-SREBP cleavage activity 1, anti-CASP3, anti-Apopain, anti-SCA-1, anti-Cysteine protease CPP32, anti-Protein Yama