

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD123/IL3RA NANOBODY®, protein (SAA1297)
PRODUCT REFERENCE	ARO-A16398
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD123/IL3RA

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-IL-3R subunit alpha, anti-IL-3 receptor subunit alpha, anti-CD123, anti-Interleukin-3 receptor subunit alpha, anti-IL3R, anti-IL-3R-alpha, anti-IL3RA, anti-IL-3RA