

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human MAGEB1 NANOBODY [®] , protein (SAA1034)
PRODUCT REFERENCE	ARO-A16400
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	MAGEB1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-MAGEXP, anti-CT3.1, anti-Cancer/testis antigen 3.1, anti-DSS-AHC critical interval MAGE superfamily 10, anti-MAGE-B1 antigen, anti-MAGE-XP antigen, anti-Melanoma-associated antigen B1, anti-DAM10, anti-MAGEB1, anti-MAGEL1