

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Anti-Human MAGEB1 NANOBODY [®] , protein (SAA1250) |
| PRODUCT REFERENCE | ARO-A16404 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | VHH-8His-Cys-tag |
| TARGET | MAGEB1 |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Homo sapiens (Human)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4) |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | anti-MAGEXP, anti-CT3.1, anti-Cancer/testis antigen 3.1, anti-DSS-AHC critical interval MAGE superfamily 10, anti-MAGE-B1 antigen, anti-MAGE-XP antigen, anti-Melanoma-associated antigen B1, anti-DAM10, anti-MAGEB1, anti-MAGEL1 |