

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human CXCL12/SDF-1 NANOBODY&reg, protein (SAA1273)
<b>PRODUCT REFERENCE</b>	ARO-A16407
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VHH-8His-Cys-tag
<b>TARGET</b>	CXCL12/SDF-1

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	anti-Interocrine reduced in hepatomas, anti-hSDF-1, anti-Stromal cell-derived factor 1, anti-SDF-1, anti-Pre-B cell growth-stimulating factor, anti-SDF1, anti-SDF1A, anti-C-X-C motif chemokine 12, anti-SDF1B, anti-CXCL12, anti-IRH, anti-PBSF, anti-hIRH