

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD154/CD40LG/TNFSF5 NANOBODY [®] protein (SAA1196)
PRODUCT REFERENCE	ARO-A16414
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD154/CD40LG/TNFSF5

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-sCD40L, anti-Tumor necrosis factor ligand superfamily member 5, anti-T-cell antigen Gp39, anti-CD40L, anti-TNFSF5, anti-CD154, anti-CD40-L, anti-TNF-related activation protein, anti-CD40 ligand, anti-CD40LG, anti-TRAP