

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human MMP12 NANOBODY [®] protein (SAA1293)
PRODUCT REFERENCE	ARO-A16415
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	MMP12

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Macrophage elastase, anti-HME, anti-ME, anti-MMP12, anti-Matrix metalloproteinase-12, anti-MME, anti-hME, anti-MMP-12, anti-Macrophage metalloelastase