

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD252/TNFSF4/OX40L NANOBODY [®] , protein (SAA1277)
PRODUCT REFERENCE	ARO-A16428
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD252/TNFSF4/OX40L

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-CD252, anti-OX40L, anti-TXGP1, anti-OX40 ligand, anti-Glycoprotein Gp34, anti-TAX transcriptionally-activated glycoprotein 1, anti-TNFSF4, anti-Tumor necrosis factor ligand superfamily member 4