

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human CD182/CXCR2 NANOBODY&reg, protein (SAA1360)
<b>PRODUCT REFERENCE</b>	ARO-A16434
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VHH-8His-Cys-tag
<b>TARGET</b>	CD182/CXCR2

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	anti-IL8RB, anti-GRO/MGSA receptor, anti-C-X-C chemokine receptor type 2, anti-CD182, anti-CXCR-2, anti-IL-8R B, anti-CXC-R2, anti-IL-8 receptor type 2, anti-High affinity interleukin-8 receptor B, anti-CDw128b, anti-CXCR2