

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CCL3/MIP-1-alpha NANOBODY [®] protein (SAA1270)
PRODUCT REFERENCE	ARO-A16439
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CCL3/MIP-1-alpha

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Tonsillar lymphocyte LD78 alpha protein, anti-CCL3, anti-G0S19-1, anti-MIP1A, anti-Macrophage inflammatory protein 1-alpha, anti-LD78-alpha(4-69), anti-SCYA3, anti-SIS-beta, anti-Small-inducible cytokine A3, anti-C-C motif chemokine 3, anti-G0/G1 swi