

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human LMO2 NANOBODY [®] , protein (SAA0919)
PRODUCT REFERENCE	ARO-A16441
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	LMO2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-RBTNL1, anti-Rhombotin-2, anti-T-cell translocation protein 2, anti-LMO-2, anti-Cysteine-rich protein TTG-2, anti-TTG2, anti-LIM domain only protein 2, anti-LMO2, anti-RHOM2, anti-RBTN2