

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD11b/ITGAM NANOBODY [®] protein (SAA1302)
PRODUCT REFERENCE	ARO-A16447
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD11b/ITGAM

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-CR-3 alpha chain, anti-Cell surface glycoprotein MAC-1 subunit alpha, anti-Leukocyte adhesion receptor MO1, anti-ITGAM, anti-CD11b, anti-CD11 antigen-like family member B, anti-Integrin alpha-M, anti-CR3A, anti-CD11B, anti-Neutrophil adherence recept