

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD326/EPCAM NANOBODY®, protein (SAA1341)
PRODUCT REFERENCE	ARO-A16451
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD326/EPCAM

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Epithelial glycoprotein 314, anti-hEGP314, anti-Major gastrointestinal tumor-associated protein GA733-2, anti-EGP314, anti-TACSTD1, anti-Adenocarcinoma-associated antigen, anti-TROP1, anti-EPCAM, anti-Epithelial glycoprotein, anti-KS 1/4 antigen, ant