

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CCL5/RANTES NANOBODY®, protein (SAA1271)
PRODUCT REFERENCE	ARO-A16456
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CCL5/RANTES

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-TCP228, anti-C-C motif chemokine 5, anti-CCL5, anti-T cell-specific protein P228, anti-D17S136E, anti-Small-inducible cytokine A5, anti-SIS-delta, anti-SCYA5, anti-T-cell-specific protein RANTES, anti-Eosinophil chemotactic cytokine, anti-EoCP