

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD9 NANOBODY [®] , protein (SAA0905)
PRODUCT REFERENCE	ARO-A16469
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD9

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Motility-related protein, anti-Leukocyte antigen MIC3, anti-TSPAN29, anti-Cell growth-inhibiting gene 2 protein, anti-CD9, anti-MIC3, anti-p24, anti-Tspan-29, anti-MRP-1, anti-5H9 antigen, anti-Tetraspanin-29, anti-CD9 antigen