

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human PSIP1 NANOBODY [®] , protein (SAA1182)
PRODUCT REFERENCE	ARO-A16482
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	PSIP1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Transcriptional coactivator p75/p52, anti-CLL-associated antigen KW-7, anti-PSIP1, anti-DFS70, anti-Dense fine speckles 70 kDa protein, anti-DFS 70, anti-PC4 and SFRS1-interacting protein, anti-LEDGF, anti-PSIP2, anti-Lens epithelium-derived growth f