

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-General MTX/Methotrexate NANOBODY [®] protein (SAA1192)
PRODUCT REFERENCE	ARO-A16489
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	MTX/Methotrexate

PRODUCT INFORMATION

ORIGIN SPECIES	<i>General</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-4-amino-10-methylfolic acid, anti-4-amino-N10-methylpteroylglutamic acid, anti-Methotrexate, anti-METHOTREXATE, anti-MTX, anti-N-[4-[[[2,4-diamino-6-pteridinyl)methyl]methylamino]benzoyl]-L-glutamic acid