

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-General AX/Amoxicillin NANOBODY [®] protein (SAA1371)
PRODUCT REFERENCE	ARO-A16492
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	AX/Amoxicillin

PRODUCT INFORMATION

ORIGIN SPECIES	<i>General</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-(2S,5R,6R)-6-[[[(2R)-2-amino-2-(4-hydroxyphenyl)acetyl]amino]-3,3-dimethyl-7-oxo-4-thia-1-azabicyclo[3.2.0]heptane-2-carboxylic acid, anti-6-(p-hydroxy- β -aminophenylacetamido)penicillanic acid, anti- β -amino-p-hydroxybenzylpenicillin, anti-Amoxicilli