

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-General AZT/Aztreonam NANOBODY®, protein (SAA1370)
PRODUCT REFERENCE	ARO-A16495
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	AZT/Aztreonam

PRODUCT INFORMATION

ORIGIN SPECIES	<i>General</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-(Z,)-2-((((2-Amino-4-thiazolyl)(((2S,3S)-2-methyl-4-oxo-1-sulfo-3-azetidiny)l)carbamoyl)methylene)amino)oxy)-2-methylpropionic acid, anti-AZT