

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD156c/ADAM10 NANOBODY [®] protein (SAA1217)
PRODUCT REFERENCE	ARO-A16501
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD156c/ADAM10

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-ADAM 10, anti-Mammalian disintegrin-metalloprotease, anti-CD156c, anti-MADM, anti-ADAM10, anti-Kuzbanian protein homolog, anti-KUZ, anti-CDw156, anti-Disintegrin and metalloproteinase domain-containing protein 10