

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human AURKA NANOBODY [®] , protein (SAA1200)
PRODUCT REFERENCE	ARO-A16504
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	AURKA

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Serine/threonine-protein kinase 15, anti-Aurora/IPL1-related kinase 1, anti-AIK, anti-ARK1, anti-ARK-1, anti-STK15, anti-Aurora kinase A, anti-STK6, anti-Aurora 2, anti-hARK1, anti-Serine/threonine-protein kinase 6, anti-AYK1, anti-AURKA, anti-Breast