

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |  |
|--------------------------|--|
| <b>PRODUCT NAME</b>      | Anti-Human BIRC5 NANOBODY <sup>®</sup> protein (SAA1237) |
| <b>PRODUCT REFERENCE</b> | ARO-A16516   |
| <b>CLONALITY</b>         | Monoclonal Antibody                                      |
| <b>ISOTYPE</b>           | VHH-8His-Cys-tag   |
| <b>TARGET</b>            | BIRC5  |

### PRODUCT INFORMATION

---

|                                |   |
|--------------------------------|---|
| <b>ORIGIN SPECIES</b>          | <i>Homo sapiens (Human)</i>   |
| <b>FORM</b>                    | Liquid  |
| <b>BUFFER</b>                  | 0.01M PBS (pH 7.4)  |
| <b>STABILITY &amp; STORAGE</b> | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| <b>ALIAS</b>                   | anti-BIRC5, anti-IAP4, anti-Apoptosis inhibitor survivin, anti-Apoptosis inhibitor 4, anti-API4, anti-Baculoviral IAP repeat-containing protein 5   |