

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-P.gingivalis fimA/Fimbrilin NANOBODY®, protein (123-123-10)
PRODUCT REFERENCE	ARO-A16526
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	P.gingivalis fimA/Fimbrilin

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Porphyromonas gingivalis</i> (strain ATCC 33277 / DSM 20709 / CIP 103683 / JCM12257 / NCTC 11834 / 2561)
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Major fimbrium subunit FimA type-1, anti-Fimbrillin, anti-Major fimbrial subunit protein type I, anti-fimA, anti-PGN_0180, anti-Bacteria